



# CLINICAL TRIAL AND REAL-WORLD EVIDENCE OF ctDNA MONITORING TO PREDICT RECURRENCE IN PATIENTS WITH RESECTED COLORECTAL CANCER

Volume 6,100076, December 2024- ESMO Real World Data and Digital Oncology ([esmorwd.org](http://esmorwd.org))

H. T. Nguyen, V.-A. Nguyen Hoang, T. V. Nguyen, N. A. L. Trinh, T. H. Pham, D. N. Le, H. H. Ho, T. D. Nguyen, H. D. Nguyen, T. S. L. Thi, N. Nguyen, D. H. Tran, M. T. Le, T. C. Dinh, T. S. Nguyen, K. C. N. The, H. Mai, M. T. Chu, D. H. Pham, N. H. T. Phuc, D. N. Vinh, D.-N. Nguyen, X.-V. Nguyen, D. S. Nguyen, T. T. D. Thi, H. Giang, H.-N. Nguyen & L. N. Tu

<sup>1</sup>University Medical Center Ho Chi Minh City, Ho Chi Minh City; <sup>2</sup>Medical Genetics Institute, Ho Chi Minh City; <sup>3</sup>Gene Solutions, Ho Chi Minh City; <sup>4</sup>Thu Duc City Hospital, Ho Chi Minh City; <sup>5</sup>Dong Nai General Hospital, Dong Nai; <sup>6</sup>Military Hospital 175, Ho Chi Minh City; <sup>7</sup>Thanh Nhan Hospital, Hanoi; <sup>8</sup>Nam Sai Gon International General Hospital, Ho Chi Minh City, Vietnam

A new study published in ESMO Real World Data and Digital Oncology demonstrated that:

- ✔ Post-operative ctDNA was a strong **independent prognostic biomarker**.
- ✔ **Performance** of ctDNA was comparable between clinical trial and real-world cohorts.
- ✔ ctDNA was **more sensitive and specific** than CEA to detect recurrence.
- ✔ High frequency of co-mutations in **TP53, APC and KRAS** was observed in relapsed patients.

## Background and Solutions for Minimal Residual Disease (MRD) Detection


### Current Limitations

- **Limited sensitivity and specificity** characterize conventional methods like ultrasound, CT scans, physical exams, and CEA blood tests, making early detection of MRD challenging
- **Lacks clinical validation and adoption**, ctDNA monitoring in cancer management is limited in developing countries.

### The ctDNA Solutions

- **ctDNA is a promising biomarker** for detecting and quantifying MRD, with strong prognostic value for cancer recurrence.
- Analysis from clinical trials and real-world data shows ctDNA as an **independent and strong prognostic biomarker for predicting CRC recurrence**, aiding personalized intervention and surveillance.


## Introducing K-TRACK™ – A Personalized ctDNA Analysis Workflow



## K-TRACK™ Performance in Clinical Trials


### 80.0% sensitivity at the post-op landmark

(2-4 weeks post-surgery) was observed. ctDNA was detected in 19.4% of patients and was associated with a significantly shorter DFS, while CEA status could not predict DFS. ctDNA detection increased the risk of recurrence by 30.3 times



### 100.0% sensitivity at post-op surveillance


was observed. ctDNA continued to show the strongest prognostic value, with a hazard ratio of 173.6 for predicting recurrence. ctDNA detection increased the risk of recurrence by 173.6 times



### 8 months Lead Time for Recurrence Detection

The study reported a median lead time of 8 months between ctDNA detection and the clinical diagnosis of recurrence.


This early detection window offers a critical opportunity for intervention before traditional imaging would detect a recurrence.



## K-TRACK™ Performance in Real-world Practice

### Infrequent Surveillance Testing:

In real-world practice, regular ctDNA surveillance was not routinely performed. Most patients were tested only once or twice post-surgery, particularly around **landmark time points** (immediately after surgery or after adjuvant therapy). This highlights potential missed opportunities for early detection of recurrence and treatment modification.




### ctDNA detection identified 87.5% of recurrences before clinical diagnosis


This early detection provides oncologists with a crucial tool for timely intervention, potentially improving patient outcomes by catching recurrence earlier than imaging methods.

### 90.0% sensitivity and 100.0% specificity for detecting recurrence

This underscores ctDNA's reliability as a Prognostic biomarker for identifying patients at high risk of recurrence. Lead time up to 5.9 months




### Suspected cases (n=8) Concordance: 87.5%



## Mutational Insights

Common mutations detected in ctDNA among relapsed patients were in the TP53, APC, and KRAS genes


These mutations, detected in both clinical trials and real-world practice, could help stratify patients for targeted therapies and further personalize treatment approaches



## Comparison to Conventional Markers

Before surgery, ctDNA detection was 91.2%, compared to 42.9% for CEA.


ctDNA detection was far more effective than carcinoembryonic antigen (CEA), the traditional blood marker for CRC.



## CASE STUDY - SILENT SIGNALS:

### ctDNA's Early Warning of Recurrence in Colorectal Cancer

#### Patient Information



#### Case ZMC111

62 years old – Female  
Sigmoid colon adenocarcinoma,  
Grade 2, stage III/IV pT4N0Mx

A 62-year-old woman with **stage III/IV sigmoid colon adenocarcinoma** underwent surgery and XELOX chemotherapy. Her **CEA and ctDNA levels normalized post-surgery**, indicating remission. However, at the **12-month follow-up, ctDNA reappeared at 0.10%, despite normal CEA and no symptoms**. A CT scan at **20 months revealed metastasis**, showing that ctDNA monitoring provided an early, non-invasive alert of recurrence **7.7 months before imaging detected it**.

#### Pre-Operative

#### Therapy Treatment

#### Monitoring

#### Recurrence

Surgery

Adjuvant  
Therapy

6.5  
months

12  
months

20  
months

CEA +  
8.1 ng/ml

Curative  
-intent

XELOX

CEA -


CEA -

Inferior  
mesenteric lymph  
node metastasis  
detected via CT  
scan

ctDNA -

ctDNA +  
0.1%

#### Recurrence Detection



Longitudinal ctDNA monitoring can serve as an early, non-invasive indicator of cancer recurrence, detecting it significantly earlier than traditional imaging methods